
# From Eye to Insight





## DMi1 – Designed for routine cell and tissue culture

#### DMi1 for quick check on cell status

For basic usage, high quality and easy-to-use product with little maintenance need.

- > S80/0.30 condenser (80 mm working distance).
- > 10x to 20x to 40x phase contrast with one hand move; automatically light intensity change between brightfield and phase contrast.
- > LED technology and a two-hour auto-off-function reduce energy costs.



Quick check on cell status

#### DMi1 with camera for cell check and documentation

For cell culture labs with limited space or high cleanness demanding and documentation required.

- > S40/0.45 condenser (40-50 mm working distance) for optimal resolution with types of dishes, flasks and microtiters with extra 10mm free space by simply moving up.
- > 12MP integrated FLEXACAM C1 camera turns your microscope into a stand-alone digital imaging station with no need for a PC.
- > Fits narrow into a high cleanness demanding cell culture room.
- > Objective guide with universal holder for an ease-of-use operation.



Cell check and documentation





Leica Microsystems CMS GmbH | Ernst-Leitz-Strasse 17-37 | D-35578 Wetzlar (Germany)

Tel. +49 (0) 6441 29-0 | F +49 (0) 6441 29-2599

www.leica-microsystems.com

### CONNECT WITH US!

