

# ROFIN FLARE+2

## LONG CONTINUOUS MAXIMUM OUTPUT

200 minutes at full power

## POWERFUL OUTPUT

most powerful optical output available

## MAINS OR CAR CHARGING

heads can be either directly attached or remote

## WATERPROOF

submersible to 300ft

## LIGHT WEIGHT

weighs only 700gms

## PURE LIGHT OUTPUT

filtered to 10<sup>-6</sup>

## SIMPLE OPERATION

turn switch handle easily turns on and off

## FULLY INTEGRATED UNIT

head and battery are fully encapsulated in one body

## 30% MORE POWERFUL

more powerful than previous generations



LEDs	7 high intensity LEDs
LEDs Power Output	up to 6500 mW
White LED Output	3300 lumens
Runtime	up to 2,000 minutes
Battery life at full power	200 minutes
Charge time	60-150 minutes

	OUTPUT	10% BAND WIDTH	PEAK
	White	400-700	-
	UV	360-385	365
	Violet	405-420	415
	Blue	435-465	450
	Cyan	485-515	505
	Green	510-545	530
	Green/Yellow	530-560	545
	Orange	590-605	595
	Red	615-635	620
	IR	835-865	850

## FLARE+ 2 SINGLE HEADS

The next generation of Forensic Light Sources: more powerful, longer battery life, smaller and lighter. The sleek design of the Flare Plus2 fits comfortably into the hand and allows the thumb to control the multi-function reed switch collar. The reed switch controls ON/OFF, intensity of output (5 steps), strobing, safety locking and SOS. The unit is coated with a fluoro-polymer for easy cleaning. It is brighter than any other hand held battery system. Lithium batteries can last up to 33 hours per charge (on 35% power) and 5 hours on full power. The UV Head produces 2450mW optical light and the 415nm head produces 6500mW. This device is waterproof and is submersible up to 300 feet.

- A-20001 UV LED 365nm
- A-20002 Purple LED 415nm
- A-20003 Blue LED 450nm
- A-20004 Cyan LED 505nm
- A-20005 Green LED 530nm
- A-20006 Green LED 545nm
- A-20024 Orange LED 595nm
- A-20007 Red LED 620nm
- A-20008 White LED 400-700nm
- A-20031 IR LED 850nm





Charge up to 4 Flashlights at a time



Pelican case w/ built in charger



## ROFIN FLARE PLUS 2 ALTERNATE LIGHT SOURCE

The next generation of Forensic Light Sources: more powerful, longer battery life, smaller and lighter. The sleek design of the Flare Plus2 fits comfortably into the hand and allows the thumb to control the multi-function reed switch collar. The reed switch controls ON/OFF, intensity of output (5 steps), strobing, safety locking and SOS. The unit is coated with a fluoro-polymer for easy cleaning. It is brighter than any other hand held battery system. Lithium batteries can last up to 33 hours per charge (on 35% power) and 5 hours on full power. The UV Head produces 2450mW optical light and the 415nm head produces 6500mW. This device is waterproof and is submersible up to 300 feet.

### KIT CONTENTS

Pelican™ case w/ built in charger  
365nm UV LED Head & Battery  
415nm Purple LED Head & Battery  
450nm Blue LED Head & Battery  
505nm Cyan LED Head & Battery  
530nm Green LED Head & Battery

545nm Green LED Head & Battery  
620nm Red LED Head & Battery  
400-700nm White LED Head & Battery  
Beam Shaping Filter  
(5°, 20°, 50° and 60° x 1°)  
Red, Yellow, UV and Orange Goggles

Auto Charger  
Power Cable  
User Manual

A-FL2KIT

