

# ROFIN FLARE+ 2

THE BRIGHTEST HAND HELD BATTERY OPERATED FORENSIC LIGHT



## **POWERFUL OUTPUT**

*most powerful optical output available*

## **MAINS OR CAR CHARGING**

*heads can be either directly attached or remote*

## **WATERPROOF**

*submersible to 300ft*

## **LIGHT WEIGHT**

*weighs only 700gms*

## **PURE LIGHT OUTPUT**

*filtered to  $10^{-6}$*

## **LONG CONTINUOUS MAXIMUM OUTPUT**

*200 minutes at full power*

## **SIMPLE OPERATION**

*turn switch handle easily turns on and off*

## **FULLY INTEGRATED UNIT**

*head and battery are fully encapsulated in one body*

## **30% MORE POWERFUL**

*more powerful than previous generations*

**DISTRIBUTED BY:**



**ARROWHEAD FORENSICS**

11030 Strang Line Road | Lenexa, Kansas 66215

PHONE 913.894.8388 FAX 913.894.8399

[www.arrowheadforensics.com](http://www.arrowheadforensics.com)



# ROFIN FLARE+ 2

KIT AND ACCESSORIES OPTIONS



CrimePak Side Pouch Available



Charge up to 4 Flashlights at a time



Pelican® case w/ built in charger

## 8 HEAD KIT CONTENTS

Pelican® case w/ built in charger  
365nm UV LED Head & Battery  
415nm Purple LED Head & Battery  
450nm Blue LED Head & Battery  
505nm Cyan LED Head & Battery  
530nm Green LED Head & Battery  
545nm Green LED Head & Battery  
620nm Red LED Head & Battery  
400-700nm White LED Head & Battery  
Beam Shaping Filter  
(5°, 20°, 50° and 60° x 1°)  
Red, Yellow, UV and Orange Goggles  
Ball Mount  
Auto Charger  
Power Cable  
User Manual

## ROFIN FLARE PLUS2 ALTERNATE LIGHT SOURCE

The next generation of Forensic Light Sources: more powerful, longer battery life, smaller and lighter. The sleek design of the Flare Plus2 fits comfortably into the hand and allows the thumb to control the multi-function reed switch collar. The reed switch controls ON/OFF, intensity of output (5 steps), strobing, safety locking and SOS. The unit is coated with a fluoro-polymer for easy cleaning. It is brighter than any other hand held battery system. Lithium batteries can last up to 33 hours per charge (on 35% power) and 5 hours on full power. The UV Head produces 2450mW optical light and the 415nm head produces 6500mW. This device is waterproof and is submersible up to 300 feet.

## ADDITIONAL KIT OPTIONS

We custom build Flare+ 2 kits to meet you specifications. Kits range from single head to 8 heads and can be packed in our CrimePak Soft Pouch cases or Pelican Watertight cases.



## ACCESSORIES

A-20028 Desk Charging System  
A-20029 Car Charging System  
A-20027 Case to Wall Charging System  
A-20011 Ball Mount  
A-20013 5° Beam Shaping Filter  
A-20014 20° Beam Shaping Filter  
A-20015 50° Beam Shaping Filter  
A-20016 60° x 1° Beam Shaping Filter