



A-FRD-001SS

- Special purpose dryer for swabs
- Horizontal airflow geometry
- Swab holder included
- 2 sizes available - 15 and 30 capacity



A-FRD-001MSS

# SAFESWAB DRYING CABINET

## EVIDENCE DRYING CABINET FOR SWABS

The HEPA main filter is formulated to contain any particulates or biological material that might be present on swabs. Safeswab is used to safely dry evidence swabs while protecting them from cross-contamination. The Safeswab Swab Dryer measures 16" W x 20" D x 16" H and features a removable rack that holds 30 Swabs. The SafeSwab Swab Dryer Mini measures 7.5"W x 8"D x 9.5"H and features a removable rack that holds 15 Swabs.



A-FRD-001SS



(side view: left)



(side view: right)



A-FRD-001MSS



(side view: left)



(side view: right)

MODEL	Racks		Swab Capacity	External (W x D x H)	Shipping (W x D x H)	WEIGHT (lbs/Kg)	
						Net	Ship

### Safeswab™ Models

A-FRD-001SS	2	30	16" x 20" x 16" 406 x 508 x 406 mm	24" x 24" x 24" 610 x 610 x 610 mm	18 / 8.2	25 / 11.4
A-FRD-001MSS	1	15	7.5" x 8" x 9.5" 191 x 204 x 242 mm	12" x 12" x 12" 305 x 305 x 305 mm	9 / 4.1	14 / 6.4



**ARROWHEAD FORENSICS**

11030 Strang Line Road | Lenexa, Kansas 66215

PHONE 913.894.8388 FAX 913.894.8399

www.arrowheadforensics.com